

ISO VANILLA 2000mg


Batch ID:	N/A	Test ID:	5188917.0029
Reported:	19-Aug-2019	Method:	TM17
Type:	Concentrate		
Test:	Pesticides		


PESTICIDE RESIDUE

Compound	Dynamic Range (ppb)	Result (ppb)	Compound	Dynamic Range (ppb)	Result (ppb)
Acephate	46 - 2151	ND*	Malathion	46 - 2151	ND*
Acetamiprid	46 - 2151	ND*	Metalaxyl	279 - 2151	ND*
Avermectin	279 - 2151	ND*	Methiocarb	46 - 2151	ND*
Azoxystrobin	46 - 2151	ND*	Methomyl	46 - 2151	ND*
Bifenazate	46 - 2151	ND*	MGK 264 1	46 - 2151	ND*
Boscalid	279 - 2151	ND*	MGK 264 2	279 - 2151	ND*
Carbaryl	46 - 2151	ND*	Myclobutanil	279 - 2151	ND*
Carbofuran	46 - 2151	ND*	Naled	279 - 2151	ND*
Chlorantraniliprole	46 - 2151	ND*	Oxamyl	46 - 2151	ND*
Chlorpyrifos	279 - 2151	ND*	Paclobutrazol	46 - 2151	ND*
Clofentezine	46 - 2151	ND*	Permethrin	279 - 2151	ND*
Diazinon	46 - 2151	ND*	Phosmet	46 - 2151	ND*
Dichlorvos	279 - 2151	ND*	Prophos	279 - 2151	ND*
Dimethoate	46 - 2151	ND*	Propoxur	279 - 2151	ND*
E-Fenpyroximate	279 - 2151	ND*	Pyridaben	279 - 2151	ND*
Etofenprox	279 - 2151	ND*	Spinosad A	46 - 2151	ND*
Etoxazole	279 - 2151	ND*	Spinosad D	279 - 2151	ND*
Fenoxycarb	46 - 2151	ND*	Spiromesifen	46 - 2151	ND*
Fipronil	279 - 2151	ND*	Spirotetramat	279 - 2151	ND*
Flonicamid	46 - 2151	ND*	Spiroxamine 1	46 - 2151	ND*
Fludioxonil	279 - 2151	ND*	Spiroxamine 2	46 - 2151	ND*
Hexythiazox	279 - 2151	ND*	Tebuconazole	46 - 2151	ND*
Imazalil	279 - 2151	ND*	Thiacloprid	46 - 2151	ND*
Imidacloprid	46 - 2151	ND*	Thiamethoxam	46 - 2151	ND*
Kresoxim-methyl	46 - 2151	ND*	Trifloxystrobin	279 - 2151	ND*

* ND = None Detected (Defined by Dynamic Range of the method)

FINAL APPROVAL


 Sam Smith
 19-Aug-2019
 3:14 PM
 PREPARED BY / DATE


 Greg Zimpfer
 19-Aug-2019
 3:47 PM
 APPROVED BY / DATE

Testing results are based solely upon the sample submitted to Botanacor Laboratories, LLC, in the condition it was received. Botanacor Laboratories, LLC warrants that all analytical work is conducted professionally in accordance with all applicable standard laboratory practices using validated methods. Data was generated using an unbroken chain of comparison to NIST traceable Reference Standards and Certified Reference Materials. This report may not be reproduced, except in full, without the written approval of Botanacor Laboratories, LLC.